

# RECOMMENDED DRILL PROTOCOL

## Straight

IMPLANT DIAMETER	Drill Diameter [mm]	Ø1.9	Ø2.0	Ø2.3	Ø2.8
		Drill Speed [RPM]	1200-1500	900-1200	800-1000

  

IMPLANT DIAMETER	Soft Bone	2mm			
	Ø3.75	Hard Bone	2/3		

	Drill to mark osteotomy site		Drill osteotomy to implant		Drill osteotomy partially according to implant		Drill through crestal bone only
---	------------------------------	---	----------------------------	---	--	---	---------------------------------

▲ The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment.

▲ The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

## DRILLING PROTOCOL LEGEND



	Drill to mark osteotomy site		Drill osteotomy to implant		Drill osteotomy partially according to implant		Drill through crestal bone only or with countersink to prepare the crest
---	------------------------------	---	----------------------------	---	--	---	--