


# ONE-PIECE SERIES | MBI™










## RECOMMENDED STRAIGHT DRILL PROTOCOL

Drill Diameter [mm]		Ø1.2	Ø1.5	Ø2.0
Drill Speed [RPM]		1200-1500	1200-1500	900-1200
IMPLANT DIAMETER	Ø2.0	Soft Bone	2/3	↓
		Hard Bone		↓
	Ø2.4	Soft Bone	↓	2/3
		Hard Bone		↓
	Ø2.9	Soft Bone		↓
		Hard Bone		↓

	Drill to mark osteotomy site		Drill osteotomy to implant		Drill osteotomy partially according to implant
---	------------------------------	---	----------------------------	---	--

The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment. The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

## MBI COMPONENTS

					
NM-T3201 Metal Housing	NM-T3202 Extra Soft Nylon Cap	NM-T3203 Soft Nylon Cap			
					
NMCD2312 Drill Ø 1.2mm	NMCD2315 Drill Ø 1.5mm	NM-X1520 Driver Ø 1.7mm	NM-X2415 Driver Ø 1.7mm	NM-T6124 Analog	NM-T4410 Transfer