

# NORIS EZ TOOTH™

A SIMPLIFIED, PREDICTABLE  
SINGLE-TOOTH IMPLANT SOLUTION.

Engineered for long-term stability  
and restorative precision.



# NORIS EZ TOOTH™ concept

Noris Single-Tooth Solution™ is a comprehensive, clinically driven approach to single-tooth implant therapy — from planning and surgical placement to the final restorative outcome.

**ONE SYSTEM.**  
**ONE PHILOSOPHY.**  
**PREDICTABLE RESULTS.**

Why a Single-Tooth Implant Is a Clinical Challenge?

A Single-Tooth case may appear simple — but clinically, it demands the highest level of precision. There is no room for error.

Every deviation directly affects function and esthetics.

High anatomical variability between patients.

A delicate balance between primary stability, bone preservation, and soft-tissue esthetics.


True success in single-tooth implant therapy depends not on the implant alone — but on a fully integrated system working in harmony.

# TWO CONNECTIONS. ONE CLINICAL PHILOSOPHY.

Noris Medical offers two implant connection platforms for Single-Tooth restoration:

Int Hex- Internal Hex Connections

CP- Conical Platform Connections



It's all about choosing the right solution for the right case

A Complete Workflow — No Compromises

## WHY NORIS EZ-TOOTH™ SOLUTIONS:



Proven clinical experience



Solutions developed by clinicians, for clinicians



Practical, clinically oriented thinking



Seamless integration of surgery, restoration, and digital workflows

# 3 A HOLISTIC SOLUTION - FULLY INTEGRATED STAGES

1

The Implant

2

The Surgical Protocol

3

The Restorative Solution

# TUFF™ UniCon

## FEATURES AND CLINICAL BENEFITS

### Advanced Conical Connection for Maximum Stability and Precision


The TUFF™ UNICON implant features a precision-engineered conical platform designed to enhance mechanical stability, optimize load distribution, and support long-term bone and soft tissue preservation.

The conical connection creates a secure, sealed interface between implant and abutment, minimizing micro-movement and enabling predictable esthetic and functional outcomes - particularly in demanding Single-Tooth restorations.



### STRUCTURAL FEATURES

- **Precision conical connection**  
for enhanced mechanical stability
- **Platform-switching design**  
supporting crestal bone preservation
- **Optimized implant geometry**  
for controlled insertion and primary stability
- **Advanced surface treatment**  
supporting osseointegration
- Compatible with fully digital and conventional restorative workflows



### CLINICAL ADVANTAGES

- **Connection Stability**  
Conical interface reduces micro-movement and enhances long-term mechanical integrity.
- **Bone Preservation**  
Platform-switching and connection design help maintain crestal bone levels.
- **Soft Tissue Control**  
Supports stable soft tissue architecture and predictable esthetic outcomes.
- **Prosthetic Precision**  
Provides accurate positioning and long-term restorative reliability.

# TUFF™ UniCon

## CP | SINGLE-TOOTH DRILLING PROTOCOLS

### RECOMMENDED DRILL PROTOCOL


#### Straight and Step Drills

Drill Diameter [mm]	Ø1.9	Ø1.5	Ø2.0	Ø2.8	Ø3.2	Ø3.65	Ø4.2	Ø5.2	ØCS 5-6
Drill Speed [RPM]	1200-1500	900-1200	900-1200	800-1000	500-700	400-700	400-600	400-600	400-600
IMPLANT DIAMETER	Ø3.25	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow]						
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		1/3 [Down Arrow]				
	Ø3.5	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow]						
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		2/3 [Down Arrow]				
	Ø3.75	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]						
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		2/3 [Down Arrow]				
	Ø4.2	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		2/3 [Down Arrow]				
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		1/3 [Down Arrow]				
	Ø5.0	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		1/3 [Down Arrow]				
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		1/3 [Down Arrow]		Ø5 [Down Arrow]		
	Ø5.5	Soft Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		2/3 [Down Arrow]		1/3 [Down Arrow]		
		Hard Bone	[Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow] → [Down Arrow]		1/2 [Down Arrow]				
		[Down Arrow]	[Down Arrow]	[Down Arrow]	2/3 [Down Arrow]	1/3 [Down Arrow]	[Down Arrow]		
		Drill to mark osteotomy site	Drill osteotomy to implant	Drill osteotomy partially according to implant	Drill through crestal bone only				

▲ The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment.

▲ The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

### Drilling Protocol Legend



# TUFF™ UniCon

## CP- SINGLE-TOOTH RESTORATION PROTOCOLS


# CONICAL CONNECTION

## SINGLE-TOOTH RESTORATION PROTOCOLS


Advanced restorative solutions for Single-Tooth implants with TUFF™ Int Hex — designed to integrate seamlessly with digital workflows, existing clinical equipment, and the level of control preferred by the clinician.

The TUFF™ UNICON conical platform offers a modern, precision-driven restorative concept for Single-Tooth implants. Designed to support both standard and highly customized restorative workflows, the system provides exceptional connection stability, soft-tissue control, and long-term clinical predictability.


With a conical connection engineered for optimal seal and load distribution, TUFF™ UNICON enables clinicians and laboratories to achieve precise, esthetic, and durable restorative outcomes across a wide range of clinical scenarios.



## Ti Base workflow



## Ti Pre-Milled Blank workflow




# TUFF™

## FEATURES AND CLINICAL BENEFITS


### TUFF™ Implant — Proven Connection. Trusted Clinical Performance.

The TUFF™ Int Hex implant system is built on a clinically validated connection concept, offering mechanical reliability, restorative flexibility, and long-term clinical success. Its Int hex geometry provides a precise and stable interface, supporting a wide range of prosthetic solutions while maintaining predictable surgical and restorative workflows.



#### STRUCTURAL FEATURES

- Proven Int hex connections design.
- Optimized implant macro-geometry for primary stability
- Precision machining for accurate implant–abutment fit
- Compatible with a wide range of prosthetic solutions
- Designed for both conventional and digital workflows



#### CLINICAL ADVANTAGES

- **Proven Reliability**  
Int hex connections with extensive clinical validation.
- **Restorative Flexibility**  
Supports multiple restorative options and workflows.
- **Clinical Predictability**  
Enables consistent and reproducible clinical outcomes.
- **Workflow Compatibility**  
Integrates seamlessly with existing clinical and laboratory systems.

### RECOMMENDED DRILL PROTOCOL

#### Straight Drill Protocol

		Drill Diameter [mm]	Ø1.9	Ø2.0	Ø2.8	Ø3.2	Ø3.65	Ø4.2	Ø5.2	ØCS 5-6						
		Drill Speed [RPM]	1200-1500	900-1200	800-1000	500-700	400-700	400-600	400-600	400-600						
IMPLANT DIAMETER	Ø3.3	Soft Bone	▼	→	▼											
		Hard Bone	▼	→	▼	→	▼	→	1/3	▼						
	Ø3.75	Soft Bone	▼	→	▼	→	▼									
		Hard Bone	▼	→	▼	→	▼	→	2/3	▼						
	Ø4.2	Soft Bone	▼	→	▼	→	▼	→	2/3	▼						
		Hard Bone	▼	→	▼	→	▼	→	▼	→	1/3	▼				
	Ø5.0	Soft Bone	▼	→	▼	→	▼	→	▼	→	1/3	▼				
		Hard Bone	▼	→	▼	→	▼	→	▼	→	1/3	▼	→	Ø5	▼	
	Ø6.0	Soft Bone	▼	→	▼	→	▼	→	▼	→	▼	→	1/3	▼		
		Hard Bone	▼	→	▼	→	▼	→	▼	→	▼	→	▼	→	▼	→

#### Step Drill Protocol

		Drill Diameter [mm]	Ø1.9	Ø2.0	Ø2.8	Ø3.2	Ø3.65	Ø4.2	Ø5.2	ØCS 5-6						
		Drill Speed [RPM]	1200-1500	900-1200	800-1000	500-700	400-700	400-600	400-600	400-600						
IMPLANT DIAMETER	Ø3.3	Soft Bone	▼	→	▼											
		Hard Bone	▼	→	▼	→	▼	→	2/3	▼						
	Ø3.75	Soft Bone	▼	→	▼	→	▼									
		Hard Bone	▼	→	▼	→	▼	→	▼							
	Ø4.2	Soft Bone	▼	→	▼	→	▼	→	▼							
		Hard Bone	▼	→	▼	→	▼	→	▼	→	2/3	▼				
	Ø5.0	Soft Bone	▼	→	▼	→	▼	→	▼	→	2/3	▼				
		Hard Bone	▼	→	▼	→	▼	→	▼	→	2/3	▼	→	Ø5	▼	
	Ø6.0	Soft Bone	▼	→	▼	→	▼	→	▼	→	▼	→	2/3	▼		
		Hard Bone	▼	→	▼	→	▼	→	▼	→	▼	→	▼	→	▼	→

▼	Drill to mark osteotomy site	▼	Drill osteotomy to implant	▼	Drill osteotomy partially according to implant	▼	Drill with countersink to prepare the crest
---	------------------------------	---	----------------------------	---	--	---	---

The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment. The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

**TUFF™**

INT HEX- SINGLE-TOOTH RESTORATION PROTOCOLS


# INT HEX CONNECTION

## SINGLE-TOOTH RESTORATION PROTOCOLS


Advanced restorative solutions for Single-Tooth implants with TUFF™ — designed to integrate seamlessly with digital workflows, existing clinical equipment, and the level of control preferred by the clinician.




## CEREC workflow



## Ti Base workflow



## Ti Pre-Milled Blank workflow



# ONE IMPLANT

MULTIPLE RESTORATIVE PATHS.  
ONE INTEGRATED SYSTEM.

Contact  
Clinical Support



Download  
Libraries



Explore  
TUFF™ System

